



Edser

Manual: Select Template




Please go to our online APP by
opening a browser and entering the
adress <http://app.edserlabs.com>



EDSER®

Please Sign In

Username 

Password 


Login

[Did you forget your login data?](#)

Enter your credentials to log in

2

Click on "Validation"



EDSER®

Online

Home

Orders

Print Prescription Form

New Order

Search Orders

Lists

Shipping

Shipment request

Chat

Orders Management

Search Orders

Orders > Search Orders

Search Orders

Order Number

Order Number

Order Entry Date

First Name (Patient)

First Name

Order Status

Select Status Order

Surname (Patient)

Surname

Product Type

Select Product Type



Second Surname (optional) (Patient)

Second Surname (optional)

Order Type

Select Order Type

Search Orders

Number	Status	Client	Patient	Material	Cover	Order Type	Entry	Exit	Price	Actions
222				3D Print 2.5mm	EVA - Marble Green	Orthotic (Modified Duplicate)				 

Showing 1 orders.

Find the order you want to edit and click on the order number to open it.


4

Click on "the red button called "Load Template", that will open a pop up window.

The screenshot shows the EDSER Orders Management interface. The top navigation bar includes 'Print', 'Resume', 'Delete', and 'Order Duplicate' buttons. The main content area is divided into several sections:

- Client Information:** Fields for Full Name, Phone, and Clinic.
- Patient Information:** Fields for Full Name, Gender, Shoe Size, Activity, Age, Weight, Shoe Width, Left Scan, and Right Scan.
- Order Information:** This section contains multiple dropdown menus and input fields for configuring the order, including:
 - Order Type: First Visit
 - Product Type: Orthotic
 - Shell Material: 3D Print 2.5mm
 - Orthotic Style: VERSATILITY
 - Top Cover: EVA - Latex 3mm (Green)
 - Bottom Cover: EVA - Black 1mm
 - Extra Cushion: Not specified
 - Length: Full Length
 - Quantity: 1
 - Medial Arch Height: Others (checked), LFT: 0 mm, RIGHT: 0 mm
 - Lateral Arch Height: Others (checked)

A red arrow points to the 'Load Template Model Data' button in the top right of the 'Order Information' section.



Now click on the arrow to open the drop down menu.

6

Select the template you want to use, that will load the data from that template. You can either use the data from the template as is or edit the order as needed until you complete the prescription.

